

# GO D PX

Digital Portable X-Ray System

rayence

rayenceusa.com



---

## Low Dose

40% Patient Dose Reduction with GreenON Wireless Panel

## Compact & Light Weight

100 lbs with DR Panel and Small Foot Print

## Affordable + Value + Reliable

Includes High Frequency Generator

## Superior Image Quality


Latest in Cesium DR Technology

Software Customized for the Mobile Application- "One Click"

## Flexibility -Convenient Tech Friendly Design

Easy to Put into Position and Move in tight spaces

---



The Rayence Portable is a powerful Medical X-Ray system that features a fully integrated Wireless DR panel and software customized for the rapid response needed in many applications while not sacrificing accuracy or image quality. Its light weight makes it easy to transport between rooms, remote sites and does not require a dedicated room normally seen with a standard Rad Room.

Max Output Power: 2.7KW

Generator Frequency: 70KHz

KV Range : 40-100kV

mA Range : 12-30mA

Input Voltage : 110/220V

Phase /Frequency : Single /50-60HZ


**GO D P X**  
options

Mobile DR Panel  
Positioning Partner



Powered by **Green ON**

Rayence innovations in manufacturing have allowed us to create a new generation of Low Dose Cesium Iodine wireless detectors to meet the needs of all applications with superior image quality at a substantially lower dose.



SNR of a Low Dose FPD is 20-25% higher than that of Normal FPD

Max 25%  
Higher SNR

Image  
Superiority

Max. 40%  
Lower Dose

Patient Dose of a Low Dose FPD is Max. 40% lower than that of Normal FPD