


rayence

DPX PLUS


Digital Mobile X-ray System



X-ray room
not required



No special
power necessary



Minimal space
requirement



Digital and
wireless



Full turnkey

rayence DPX PLUS





No X-ray Room Needed

DPX PLUS Mobile X-Ray System eliminates the need for a dedicated exam room, freeing up valuable space in your facility for other essential functions.

Standard Power Compatibility

The DPX PLUS operates with standard 120V power, eliminating the necessity for any special power arrangements. It's seamlessly compatible with your existing infrastructure.



Minimal Space Requirement

DPX PLUS' compact design, requiring just 8' x 10' of space, optimizes space utilization. It's perfect for facilities with limited space.

Digital and Wireless Technology

Experience the convenience and efficiency as images are captured digitally and transmitted wirelessly for instant access and analysis.



Full Turnkey Solution

DPX PLUS is a complete "plug and play" solution designed for the urgent care setting. Enjoy a hassle-free setup and operation, allowing you to focus on what matters most – providing exceptional care to your patients.



1417WCE-GF

The top of the line glass free model with 100µm pixel size



100µm



Accurate AED



MultiPixel



16 hours



Light Weight



IP67
Water & Dust
Resistant



Slip &
Impact
Protection

Specification

Parameter	DPX PLUS
Generator Type	High frequency Resonant Inverter Microprocessor controlled
KVP Range	40 – 120
MA Range	15 – 60 Selectable
MAS Range	0.125 – 200
Exposure Time	0.01 – 3.33 Seconds
X-ray Tube	Stationary anode 15° Target 1.2 fs 30,000 Heat Unit Storage
Collimator	200° Rotation – High output LED lamp Laser Indicators
Exposure Switch	2 Position Prep/X-ray
Mechanical Stand	Welded Aluminum construction

* Specifications subject to change without prior notice *

Dealer Information: