


rayence

DPX

Digital Mobile X-ray System



Compact &
Lightweight



No special
power necessary



X-ray room
not required



Digital and
wireless



Full turnkey

rayence DPX




Compact & Lightweight

The DPX was designed for maneuverability. Its compact and lightweight design allows the operator to effortlessly move the system wherever needed.

Standard Power Compatibility

The DPX operates with standard 120V power, eliminating the necessity for any special power arrangements. It's seamlessly compatible with your existing infrastructure.



1417WCE-HR

The top of the line high resolution model with 100 μ m pixel size



100 μ m



Accurate AED



MultiPixel



16 hours



Light Weight



IP67
Water & Dust
Resistant



Slip &
Impact
Protection

No X-ray Room Needed

DPX Mobile X-Ray System eliminates the need for a dedicated exam room, freeing up valuable space in your facility for other essential functions.



Digital and Wireless Technology

Experience the convenience and efficiency as images are captured digitally and transmitted wirelessly for instant access and analysis.



Specification

Parameter	DPX
Generator Type	High frequency Resonant Inverter Microprocessor controlled
KVP Range	40 – 100
MA Range	15 & 30 Selectable
MAS Range	0.125 – 120
Exposure Time	96 steps
X-ray Tube	Stationary anode 15° Target 1.0 fs 30,000 Heat Unit Storage
Collimator	200° Rotation – High output LED lamp Laser Indicators
Exposure Switch	2 Position Prep/X-ray 15' Cable
Mechanical Stand	Welded Aluminum construction

* Specifications subject to change without prior notice *

Dealer Information:

Full Turnkey Solution

DPX is a complete "plug and play" solution designed for the urgent care setting. Enjoy a hassle-free setup and operation, allowing you to focus on what matters most – providing exceptional care to your patients.

